
PolySNAP M Version 2.0.1

- Fixed an issue in countries using the comma as a decimal point where shifted patterns were not being displayed correctly in the Numerical Results pane.
- A new program splash screen and window background are displayed.
- The program version number is now displayed in program title bar.
- The program, documentation and website now use the new Bruker logo.
- A new option has been added to force quantitative analysis to be attempted on all input samples regardless of how similar they are to the reference phases. This can be selected via the '*Perform analysis on all samples regardless of similarity*' checkbox in the Options: Automatic: Advanced / Quantitative options dialog box.

PolySNAP M Version 2.0.0

- The installer now adds a copy of the tutorial files to a folder named *tutorial* in the SNAP Program Files folder.
- The copy-protection system has been updated, and the *License* menu item removed as it is no longer used.
- PolySNAP M incorporates a limited subset of the PolySNAP style graphical analysis and cluster analysis tools, for datasets of up to 96 patterns.
- The program manual and other documentation has been comprehensively reorganised.
